

Anterior Watson Jones Total Hip Arthroplasty System

Instrument system specifically designed
for Direct Anterior approach THR

PRODUCT NO.'S:

6300-00 [Complete Set]

Included in Set or available Individually:

6301-L [Anterior Watson Jones Awl – Left]

1 Overall Length: 12.375" (31,5 cm)

6301-R [Anterior Watson Jones Awl – Right]

2 Overall Length: 12.375" (31,5 cm)

6302-01 [Anterior Watson Jones Lighted Mueller Retractor]

3 Blade Width: 27 mm

Handle Length: 5.5" (14 cm)

Overall Length: 12.75" (32,4 cm)

6303-01 [Anterior Watson Jones Wide Lighted Retractor]

4 Blade Width: 40 mm

Blade Tip Length: 17 mm

Overall Length: 12.75" (32,4 cm)

6304-01 [Anterior Watson Jones Narrow Lighted Retractor]

5 Blade Width: 26 mm

Blade Tip Length: 23 mm

Overall Length: 16.5" (41,9 cm)

6305 [Anterior Watson Jones 90° Cobra Retractor]

6 Blade at Widest: 37 mm

Blade Depth: 7.5" (19,1 cm)

Blade Prong: 18 mm Long X 10 mm Wide

Overall Length: 12.5" (31,8 cm)

Handle Length: 9.5" (24,1 cm)

6306 [Anterior Watson Jones Deep Hohmann Retractor]

7 Blade Width: 17 mm

Prong Length: 34 mm

Overall Length: 9.25" (23,5 cm)

6307 [Anterior Watson Jones Straight Hohmann Retractor]

8 Prong Length: 4,2 cm

Overall Length: 9.5" (24,1 cm)

6308 [Anterior Watson Jones Femoral Starter Drill]

9 Drill Tip Dimensions: 18 mm Long X 12 mm Diameter

Overall Length: 6" (15,2 cm)

Shaft Length: 5.25" (13,3 cm)

Optional Instrument (Not Included In Set):

6303 [Anterior Watson Jones Wide Retractor without Light]

Blade Width: 40 mm

Blade Tip Length: 17 mm

Overall Length: 12.75" (32,4 cm)

Lighted retractors attach to a fiber optic light cable with
ACMI (female) connector and can be steam sterilized.



Optional quick-connect driver (NOT INCLUDED
IN RETRACTOR SET) for use with the starter drill:

PRODUCT NO.:


8248 [Fixed Driver with
Zimmer Hall Quick-connect]

Overall Length: 5.75" (15,6 cm)

Handle Width: 4.625" (11,6 cm)



Fixed Driver
with Zimmer Hall
Quick-connect



1 Awl – Left
#6301-L

2 Awl – Right
#6301-R

3 Lighted Mueller Retractor
#6302-01

4 Wide Lighted Retractor
#6303-01

5 Narrow Lighted Retractor
#6304-01

6 90° Cobra Retractor
#6305

7 Deep Hohmann Retractor
#6306

8 Straight Hohmann Retractor
#6307

9 Femoral Starter Drill with
Zimmer Hall Quick-connect
#6308



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMEDI



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362