

# Argintar Claw Drill Guide Wire/Suture Passer

Designed by Evan Argintar MD

*Expandable claw design allows for minimally invasive, reproducible one-step wire/suture passage*

Especially helpful during quad/patella tendon repair surgery and other applications where a suture will be passed—particularly when soft tissue dissection is to be minimized, such as multi-ligament knee reconstruction (MCL/LCL), elbow reconstruction (ULCL/MCL), foot-ankle reconstruction (ATFL), and wrist reconstruction (DRUJ).

**PRODUCT NO:**

8315-00 [Set: (1) Claw, (1) Wire/Suture Pin]


8315-01 [Claw Unit]

Maximum Internal Opening: 2.5" (6,4 cm)

Product Dimensions: 2.5" x 4" (6,4 cm x 10,2 cm)

1227 [3/32" (2mm) Pin with Wire/Suture Hole]

Overall Length: 6" (15,2cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362