



# Arthroscopic Shoulder Rasp

*Used to abrade the anterior scapular neck to stimulate a vascular healing response*

Terminal bend matches the angle of the scapular neck from a standard anterior portal.

**PRODUCT NO:**

2310

Overall Length: 9.625"

Handle Length: 3.5"



**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

**Innomed-Europe** Tel. +41 41 740 67 73  
Fax +41 41 740 67 74



**1.800.548.2362**