

# Bacastow Arthroscopic Deltoid Lift Retractor

Designed by David Bacastow, MD


*Designed to help improve visualization of the subacromial space during rotator cuff repair or balloon spacer placement*



**PRODUCT NO:**

5081

Overall Length: 7.5" (19,1 cm)  
Depth from Bend: 2.25" (5,7 cm)  
Shaft Diameter: 4,7 mm  
Diameter at End: 6 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe**

Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

# INNOMED