

Bacastow Femoral Cement Osteotome


Designed by David Bacastow, MD

Uniquely shaped osteotome designed to help trim away cement from around a femoral knee component

Ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

PRODUCT NO:

5234
Overall Length: 9.25" (23,5 cm)
Width: 6,5 mm
Tongue Length: 7 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362