


Bach Graft Pusher Set

Designed by Bernard Bach Jr., MD



Double Pointed
Graft Pusher

Cupped
Graft Pusher

Bach Double Pointed Graft Pusher

*Flat head 3/8" wide by 1/8" thick
with two tines designed to help
stabilize graft when sliding
up the tibial tunnel*

Bach Double Pointed Graft Pusher
#5080-01

Bach Cupped Graft Pusher
#5080-02

Bach Cupped Graft Pusher

*Designed for femoral bone
plug placement using the
"push in" technique*

PRODUCT NO'S:	
5080-00	[Set]
Set Includes / Available Individually:	
5080-01	[Bach Double Pointed Graft Pusher]
Overall Length: 10" (25,4 cm)	
Handle Length: 4.25" (10.8 cm)	
Pronged End Dimensions: .375" x .125" (7,9 x 3,2 mm)	
Concave End Inner Depth: 2 mm	
5080-02	[Bach Cupped Graft Pusher]
Overall Length: 10" (25,4 cm)	
Handle Length: 4.25" (10.8 cm)	
Concave End Diameter: 5 mm	
Concave End Inner Depth: 2 mm	



www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362