


Beard IM Nail Guide Wire Clamp


Designed by David Beard, MD

Designed to help provide quick grasp-and-release of an IM guide wire for positioning and advancement along the length of the guide wire

- ▶ Anatomic pistol grip for comfortable use
- ▶ Facilitates fracture reduction in appropriate cases
- ▶ Universal to all fixation systems using IM guide wires and IM fixation
- ▶ For use with pins up to 4 mm



Available with or without ratchet



PRODUCT NO'S:

3019 [Beard Clamp with Ratchet]
Dimensions: 5.5" w x 6" h (14 cm x 15,2 cm)

3019-01 [Beard Clamp without Ratchet]
Dimensions: 5.5" w x 6" h (14 cm x 15,2 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362