



Beicker Hammerhead Rongeur

Designed by Clint Beicker, MD

Designed to help remove osteophytes from around the acetabulum, tibia, and glenoid



PRODUCT NO:
1775-05
Overall Length: 8" (20,3 cm)
Jaw Bite: 15 mm x 7 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

INNOMED