


Bent Hohmann Retractors—Narrow

Helps retract tissues at the margins of the joint

Useful for retracting tissues at the margins of the joint. Can be passed over the margins of the joint and held in place with weights or by hand.

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.



PRODUCT NO'S:

7110 [Narrow] Overall Length: 9.75" (23.8cm) Handle Length: 7" (17.8cm) Blade Width: 19mm Depth from Bend: 4.75" (12.1cm)
7110-R* [OrthoLucent™ Narrow] Overall Length: 9.75" (23.8cm) Handle Length: 7" (17.8cm) Blade Width: 19mm Depth from Bend: 4.75" (12.1cm)
7110-01 [Narrow with Extra Long Handle] Overall Length: 11.5" (29.2cm) Handle Length: 10" (25.4cm) Blade Width: 19mm Depth from Bend: 4.75" (12.1cm)
7111 [Narrow with Extra Grip Tip] Overall Length: 9.75" Blade Width: 19 mm Depth from Bend: 4.25"
7115 [Short-tipped Narrow] Overall Length: 8.625" (21.9cm) Handle Length: 7" (17.8cm) Blade Width: 19mm Depth from Bend: 4.4" (11.2cm)
7115-01 [Short-tipped Narrow with Extra Long Handle] Overall Length: 11" (27.9cm) Handle Length: 10" (25.4cm) Blade Width: 19mm Depth from Bend: 4.25" (10.8cm)
7115-03 [Extra Deep Short-tipped Narrow] Overall Length: 12.125" (31.1cm) Handle Length: 9.75" (24.8cm) Depth from Bend: 6.25" (15.9cm) Blade Width: 19mm

PROUDLY MADE IN THE USA * **MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND**


Short-tipped designed by Carl DiRaimondo, MD

Extra Grip Tip design modification by Alfred A. Durham, MD



Extra Grip Tip

Helps to stabilize the retractor around the bone



Narrow

OrthoLucent™ Narrow

Narrow with Extra Long Handle

New!
Narrow with Extra Grip Tip #7111

Short-tipped Narrow

Short-tipped Narrow with Extra Long Handle

Extra Deep Short-tipped Narrow

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362