


Bhargava Modular Offset Cup Liner Impactor

Designed by Tarun Bhargava, MD

Designed to help impact an acetabular cup liner during minimally invasive direct anterior and MIS posterior approach THR

- ▶ Used in conjunction with individual interchangeable heads (sold separately) which fit securely onto the impactor end
- ▶ Helps avoid edge loading and improper seating of the liner that can occur with a straight impactor
- ▶ Uses the same heads as the Innomed CupX Acetabular Cup Extraction System



Interchangeable Head(s) Sold Separately

PRODUCT NO:
5031
Overall Length: 15.85" (40.2 cm)
Platform End Diameter: 1" (2.54 cm)



INDIVIDUAL INTERCHANGEABLE STEEL HEADS	
5202-22	22mm
5202-26	26mm
5202-28	28mm
5202-32	32mm
5202-36	36mm
Optional Size:	
5202-38	38mm

