


# Bhargava Modular Offset Cup Liner Impactor

Designed by Tarun Bhargava, MD

*Designed to help impact an acetabular cup liner during minimally invasive direct anterior and MIS posterior approach THR*

- ▶ Used in conjunction with individual interchangeable heads (sold separately) which fit securely onto the impactor end
- ▶ Helps avoid edge loading and improper seating of the liner that can occur with a straight impactor
- ▶ Uses the same heads as the Innomed CupX Acetabular Cup Extraction System



Interchangeable Head(s) Sold Separately

**PRODUCT NO:**  
5031  
Overall Length: 15.85" (40.2 cm)  
Platform End Diameter: 1" (2.54 cm)



INDIVIDUAL INTERCHANGEABLE STEEL HEADS	
5202-22	22mm
5202-26	26mm
5202-28	28mm
5202-32	32mm
5202-36	36mm
<b>Optional Size:</b>	
5202-38	38mm



INDIVIDUAL INTERCHANGEABLE DELRIN HEADS	
5202-00	Set w/Case
5202-39	39mm
5202-40	40mm
5202-41	41mm
5202-42	42mm
5202-43	43mm
5202-44	44mm
5202-45	45mm
5202-46	46mm
5202-47	47mm
5202-48	48mm
5202-49	49mm
5202-50	50mm
5202-51	51mm
5202-52	52mm
5202-53	53mm
5202-54	54mm
5202-55	55mm
5202-56	56mm
5202-57	57mm
5202-58	58mm
5202-59	59mm
5202-60	60mm



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.