

Blount Knee Retractor

Designed by James B. Stiehl, MD



Designed for retraction in total knee arthroplasty

The long narrow blade easily fits above the capsular ligament at the joint line.
Can also be used for knee revision, fitting easily around the implant.

PRODUCT NO:

4850

Overall Length: 8.5" (21,6 cm)
Prong Width: 9 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362