

# Bolanos Modified Chandler Retractor

Designed by Alberto Bolanos, MD

*Used for retracting tissue away from the bone*



**PRODUCT NO:**

3222

Overall Length: 7.5" (19,1 cm)

Blade Width at Widest: 1" (2,54 cm)



USA MADE

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.