

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Universal Bone Grafting/ Impacting Forceps

The forceps are designed with grasping ends for delivery of bone graft. When the graft is in place, the forceps are closed, which forms the ends into an impacting punch. A striking platform is attached to the end of the forceps for tapping and tamping the graft. Four end diameters are available in two lengths.

Product No's:
5010-01 [6" Length, 1/8" Diameter End]
5010-02 [6" Length, 3/16" Diameter End]
5010-03 [6" Length, 1/4" Diameter End]
5010-04 [6" Length, 5/16" Diameter End]
5050-01 [10" Length, 1/8" Diameter End]
5050-02 [10" Length, 3/16" Diameter End]
5050-03 [10" Length, 1/4" Diameter End]
5050-04 [10" Length, 5/16" Diameter End]

Designed by
J. A. Amis, MD

*Bone graft can be grasped, placed & impacted
without changing hands or instruments*

