

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Multi-Size Bone Hooks

Designed for proximal femoral elevation in total hip replacement or in other surgery with a similar need for bone manipulation. The instrument has a blunt tip and a large handle to accommodate the use of two hands if desired.

Designed by
R.L. Wixson, M.D.



Product No's: 5910 [Small]
5915 [Medium]
5920 [Large]

Bone Hook with Cable/Wire Hole

Designed by
R.L. Wixson, M.D. &
J. McCarthy, M.D.

Product No: 5920-01 [Large with Cable/Wire Hole]

Used for Proximal Femoral Elevation



Ends drawn at actual size