

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Boynton Punch

The Boynton Punch is helpful in removing trial, femoral and revision total knee components. The flange end fits onto the flange of a femoral knee component or trial.

Product No's: 5120-01 [Standard]
5120-02 [Offset]

Designed by
L. Boynton, MD

Helpful in removing components

