



# Burkhead Glenoid Retractor



Designed by Wayne Burkhead, Jr, MD

*The retractor bar presses against the glenoid while the end of the retractor puts pressure on the posterior capsule*

## PRODUCT NO'S:

5839 [Large]  
Overall Length: 9.125" (23,2 cm)  
Blade Width at End: 1.5" (3,8 cm)

5839-SM [Small]  
Overall Length: 8.75" (22,2 cm)  
Blade Width at End: 1" (2,54 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

# INNOMED