



Calcaneal Spreader

Designed by Michael Forness, DO

Separates the calcaneal osteotomized bone for placement of tricortical bone graft

Pads have a large surface area, which easily separates the calcaneal osteotomized bone for placement of tricortical bone graft. Large pad surface area helps prevent the compression of soft calcaneal cancellous bone.



Smooth Pads



Grooved Pads



PRODUCT NO'S:

1880 [Smooth Pads]
Overall Length: 7" (17,8 cm)
Pad Dimensions: 15 mm x 12 mm

1881 [Grooved Pads]
Overall Length: 7" (17,8 cm)
Pad Dimensions: 15 mm x 12 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.