


Calibrated Ortho Spreader with Slotted Tips

Designed by Jason Bariteau, MD

*A lamina spreader with a very thin closed profile,
designed to enable distraction in tight joint spaces like
the subtalar and talonavicular joints*



PRODUCT NO:

1841
Overall Length: 6.75" (17,1 cm)
Prong Length: .5" (12,7 mm)
Calibrations: 10 mm to 35 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED