



Cannulated Fracture Awl

*Helps to reduce fractures without slipping off the bone,
and cannulated to allow the placement of k-wires*



Cannulated



PRODUCT NO:

8091

Overall Length: 7" (17,8 cm)

Handle Length: 3.3" (8,4 cm)

Cannula fits wire up to: .062" (1.6 mm)



USA MADE

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED