



Carpal Tunnel Release Guide and Blade Set

Guide designed by Peter J. Evans, MD, PhD

Guide designed to help protect the median nerve while providing a track that allows for the smooth advance of the blade to divide the transverse carpal ligament during a mini-open, non-endoscopic approach




Carpal Tunnel Release Blade
#1124-02 (Pack of 2)

Evans Carpal Tunnel Guide
#1128

Surgical technique available.

PRODUCT NO'S:
1124-00 [Set of One Guide and One Blade]
Set Includes:
[Carpal Tunnel Release Blade] Overall Length: 4.75" (12 cm) Blade Length: 1.375" (3.5 cm) Blade Height: .22" (5.6 mm) Blade Width: 1 mm
1128 [Evans Carpal Tunnel Guide] Overall Length: 8" (20.3 cm) Blade Guide Widths: 2 mm and 5 mm
Replacement Parts:
1124-02 [Carpal Tunnel Release Blades] Pack of 2 <i>One blade included in set, two with this number</i>
1128 [Evans Carpal Tunnel Guide]



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.