

Chandran Anterior Retractor for THR


Designed by Rama E. Chandran, MD

Design helps to expose the anterior rim of the acetabulum and helps prevent displacement of the retractor while reaming the acetabulum during direct anterior hip replacement



SAVANNAH (2020) 12345678901011
32 1161 63111111

PRODUCT NO:
63111
Overall Length: 13.5" (34,3 cm)
Depth from Flat: 4.75" (12,1 cm)
Blade Width: .625" (1,6 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



www.innomed.net

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.