


Chandran Cannulated Dilator/Sizer for Reconstruction

Designed by Rama E. Chandran, MD

Designed for dilating and sizing the bony tunnel during ACL reconstruction

Cannulated for Wires Up To .062" (1,6 mm).

Can also be used for sizing the tenodesis screws.



PRODUCT NO:
8204
Overall Length: 9" (22,9 cm)
Handle Width: 3.5" (8,9 cm)
Cannulated for Wires Up To: .062" (1,6 mm)
Sizer Markings: From 4 to 9 mm



ISO 13485:2016

www.innomed.net



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362