



Chandran Double Ball Spike

Designed by Rama E. Chandran, MD

Designed to help rotate and control a butterfly bone fragment for fixation



PRODUCT NO:
8027
Overall Length: 12.5" (31.8 cm)
Handle Length: 4.625" (1,17 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED [®]