



Chandran Modified Knee Retractor

Designed by Rama E. Chandran, MD

Designed to help protect the patellar tendon during robotic assisted total knee replacement

Also useful to retract structures on the lateral side of the tibia.



PRODUCT NO:
7117
Overall Length: 8.75" (22,2 cm)
Blade Width: .625" (1,6 cm)
Blade Depth from Handle: 5" (12,7 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED