



Charnley/Sorrells Low-Profile Hip Arthroplasty Retractor System

Designed by R. Barry Sorrells, MD

Conforms to the thigh, providing low-profile self-retaining exposure of the femur and acetabulum in Total Hip Arthroplasty

- ▶ Provides a varied selection of blades, both pre-shaped and malleable, for both superficial and deep retraction
- ▶ Curved blades conform with the thigh
- ▶ Self-retaining and easily adjustable



PRODUCT NO'S:	
7318-01 [Frame] Length: 14.5" Maximum Width: 12"	7318-04 [Medium Narrow Blade] Blade Width: 1.25" Blade Depth: 2.25" Overall Length: 8.5"
7318-02 [Small Narrow Blade] Blade Width: 1.25" Blade Depth: 1.25" Overall Length: 8.5"	7318-06 [Medium Malleable Blade] Blade Width: 1.25" Blade Depth: 2.25" Overall Length: 8.5"
7318-03 [Small Wide Blade] Blade Width: 2" Blade Depth: 1.25" Overall Length: 8.5"	7318-07 [Large Narrow Blade] Blade Width: 1.25" Blade Depth: 3.25" Overall Length: 8.5"



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362