




Chung Weitlaner Retractor

Designed by Raymond Chung, MD


Features long prongs allowing use of a small retractor in a small, but deep wound

Particularly useful for volar approach to the distal radius in large patients and for mini-open rotator cuff surgery.



3x2 Prongs

Prong lengths of 25mm and 30mm available with either sharp or blunt tips



3x4 Prongs

PRODUCT NO'S:	
Blunt Tips	Sharp Tips
5065 [2x3 Prongs] Blade Depth: 25 mm Overall Length: 4.5" (11,4 cm)	5066 [2x3 Prongs] Blade Depth: 25 mm Overall Length: 4.5" (11,4 cm)
5065-01 [3x4 Prongs] Blade Depth: 25 mm Overall Length: 4.5" (11,4 cm)	5066-01 [3x4 Prongs] Blade Depth: 25 mm Overall Length: 4.5" (11,4 cm)
5067 [2x3 Prongs] Blade Depth: 30 mm Overall Length: 4.5" (11,4 cm)	5068 [2x3 Prongs] Blade Depth: 30 mm Overall Length: 4.5" (11,4 cm)
5067-01 [3x4 Prongs] Blade Depth: 30 mm Overall Length: 4.5" (11,4 cm)	5068-01 [3x4 Prongs] Blade Depth: 30 mm Overall Length: 4.5" (11,4 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362