



Coated Inserter for Reverse Shoulder Glenosphere Components

Designed by Michael Radon, Ilya Voloshin, MD, and Nathan Mineo

Designed to aid in the insertion of glenospheres in limited exposure patients, allowing for insertion from the side, and coated to help protect marring component surfaces



PRODUCT NO:
5071
Overall Length: 9.5" (24,1 cm)
Inserter Arm Angle: 30°



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362