

Cobb Elevators

Ultra hard titanium nitride coating helps to extend blade life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.



Two Sizes Available
With or Without Teeth

New!
STYLE WITH TEETH

PRODUCT NO'S:	
WITH TEETH	WITHOUT TEETH
3432 [1/2" with Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1/2" (13 mm)	3436 [1/2" without Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1/2" (13 mm)
3434 [1" with Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1" (25,4 mm)	3438 [1" without Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1" (25,4 mm)

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362