


# Narrow Cobra Retractors

*A general purpose instrument for use around the femur and acetabulum in minimally invasive surgery*

Lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector. Retractor can be steam sterilized.



*Extra Long Narrow*

*Narrow*

*Lighted Narrow*



## Light Source Cable Adapters

PRODUCT NO'S:
8009-S [ACMI to Storz Adapter]
8009-W [ACMI to Wolf Adapter]



PRODUCT NO'S:
6120-04 [XL Narrow] Overall Length: 15.5" (39,4 cm) Handle Length: 11" (27,9 cm) Blade Width: 19 mm
6120 [Narrow] Overall Length: 11.75" (29,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width: 19 mm
6120-L [Lighted Narrow] Overall Length: 11.75" (29,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width: 19 mm



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016  
© 2023 Innomed, Inc.

103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362