



Curved Anterior Retractors



Narrow #7804-01



Narrow
#7804-01

Standard
#7804-02

Wide
#7804-03

PRODUCT NO'S:	
7804-01 [Narrow]	Overall Length: 10" (25,4 cm) Blade Depth: 5" (12,7 cm) Upper Blade Width: .6" (15 mm) Blade Width at Tip: .14" (3,5 mm)
7804-02 [Standard]	Overall Length: 10" (25,4 cm) Blade Depth: 5" (12,7 cm) Upper Blade Width: .8" (20 mm) Blade Width at Tip: .14" (3,5 mm)
7804-03 [Wide]	Overall Length: 10" (25,4 cm) Blade Depth: 5" (12,7 cm) Upper Blade Width: 1.44" (36 mm) Blade Width at Tip: .14" (3,5 mm)



Retractor Profile



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362