

DAA Canal Finder Rasp

Designed to help begin preparation of the femoral canal prior to stem broaching – features a large handle with a striking plate end



PRODUCT NO:

C1026

Overall Length: 13.125" (33,4 cm)

Handle Length: 6" (15,2 cm)

Maximum Shaft Diameter: 13 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED