



DMP Wire Tightener

Used to hand tighten a cerclage wire around a bone

Now with four wire holes – two for up to 20 gauge wires, and two for up to 18 gauge wires. T-Handle end is used to hand tighten a wire.



Designed by DMP

PRODUCT NO:

8729 [T-Handle]

Overall Length: 4.5" (11,4 cm)

Handle Width: 2.625" (6,7 cm)

End Diameter: 15 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.