

Das Anterior Hip Bolster Assembly

Design modification by Amal Das, MD of original design by Benjamin M. Frye, MD

Designed to help provide counter resistance on the contralateral hip during reaming and implant insertion in direct anterior hip arthroplasty



PRODUCT NO'S:

4166-00 [Das Anterior Hip Bolster Assembly]

Set Includes / Available Individually:

4166-01 [Das Anterior Hip Bolster Support]

4166-02 [Das Anterior Hip Bolster Rod]

Set Includes / Replacement Parts:

2595 [Table Clamp]

4150-PD2 [Positioning Pads - Set of 2]

4150-PS [Post Screw]



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362