

# Modified Double Prong Acetabular Retractors

*Retracts the femur anteriorly during total hip arthroplasty*


Designed to retract the femur anteriorly during total hip arthroplasty. It is hooked over the anterior pelvic brim. Weights can be added to assist in exposure and to help hold the retractor in place.



Standard  
#6170

Narrow  
#6175

PRODUCT NO'S:	
6170 [Standard]	Blade Width: 44 mm
	Overall Length: 12.5" (31,8 cm)
6175 [Narrow]	Blade Width: 32 mm
	Overall Length: 12.5" (31,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362