

# Deep Cobra Retractor

*A general purpose instrument for use around the femur and acetabulum in larger patients*

Lighted retractor can be attached to a fiber optic light cable with ACMI (female) connector. Retractor can be steam sterilized.



## Light Source Cable Adapters

### PRODUCT NO'S:

8009-S [ACMI to Storz Adapter]  
8009-W [ACMI to Wolf Adapter]



### PRODUCT NO'S:

6135 [Deep]  
Overall Length: 16" (40,6 cm)  
Handle Length: 7" (17,8 cm)  
Blade Width at Widest: 33 mm

6135-L [Lighted Deep]  
Overall Length: 16" (40,6 cm)  
Handle Length: 7" (17,8 cm)  
Blade Width at Widest: 33 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362