


Deep Meyerding Retractor with Ergonomic Handle

*A self-retaining soft tissue retractor for
use in hip, knee, and shoulder surgery*



PRODUCT NO:
6244
Overall Length: 8.5" (21,6 cm)
Blade Depth: 3.5" (8,9 cm)
Blade Width: 1" (2,54 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362