

Dennis Single Prong Hip Retractor

Designed by Douglas A. Dennis, MD

An acetabular retractor designed for total hip replacement



PRODUCT NO:

6028

Overall Length: 13.25" (33,7 cm)
Upright Handle Length: 6" (15,2 cm)
Blade Depth: 5" (12,7 cm)
Blade Width: 34 mm
Prong Depth: 1" (2,54 cm)
Prong Width: 4,5 mm



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362