



# Dingo Modified Humeral Head Retractor

*Designed to help lever and displace the proximal humerus posteriorly*



**PRODUCT NO:**  
**2438**  
Overall Length: 10.1" (25,7 cm)  
Blade Width: 1" (2,54 cm)  
Blade Curl Up at End: .56" (1,4 cm)



ISO 13485:2016

[www.innomed.net](http://www.innomed.net)



# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362