


Distal Femur Distractor

Helps distract the distal femur away from the proximal tibia

Designed to distract the distal femur away from the proximal tibia during total knee surgery. It is inserted into a pre-drilled hole in the distal femur. The bent handle allows the femur to be distracted away from the tibia. The intramedullary rod portion is fluted.



PRODUCT NO'S:	
4220-00 [Standard Handle]	
Overall Length: 12.75" (32,4 cm)	
Rod Offset from Handle: 4.5" (11,4 cm)	
4220-01 [Upward Bent Handle]	
Overall Length: 17.5" (49,6 cm)	
Rod Length from Bend: 12.75" (32,4 cm)	
Rod Offset from Handle: 4.5" (11,4 cm)	



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.