

Dozier Radiolucent Bennett Hip Fracture Retractor

Designed by John K. Dozier, MD

Can be kept in place while using image intensification or taking an x-ray

Designed to be used in hip fractures with the advantage that the retractor can be kept in place while using image intensification or taking an x-ray. The handle can be rotated to the right or left for surgeon preference. May be steam or gas sterilized.



Handle can be rotated to right or left

PRODUCT NO:

6870

Handle Length: 6.75" (17,1 cm)

Blade Length: 8.5" (21,6 cm)

Blade Width at Widest: 67 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362