


# Durham Offset Zelpi Retractor

Designed by Alfred Durham, MD

*Staggered depth retractor designed for exposure during total hip and total shoulder surgery*

- ▶ In hip surgery, with the handle towards the surgeon, the longer leg is on the inside.
- ▶ In shoulder surgery, with the handle downward, the longer leg is on the outside.
- ▶ The longer leg extends 1.1" (2,8 cm) deeper.



PRODUCT NO'S:
1573-L [Left] Overall Length: 8.5" (21,6 cm) Leg Depths: 3.1" & 4.2" (7,9 cm & 10,7 cm)
1573-R [Right] Overall Length: 8.5" (21,6 cm) Leg Depths: 3.1" & 4.2" (7,9 cm & 10,7 cm)

**PROUDLY  
MADE  
IN THE  
USA**

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Scan to  
Launch Our  
Website



**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**