

# Eickmann Knee Revision Set

Designed by Thomas Eickmann, MD

*Used for total knee revision*



## PRODUCT NO'S:

5470-00 [Complete Set]

### Individual Instruments:

5470-08 [8 mm Chisel]

Osteotome Width: 8 mm  
Blade Length: 2.375" (6 cm)  
Overall Length: 7.375" (18,7 cm)

5470-11 [11 mm Chisel]

Osteotome Width: 11 mm  
Blade Length: 2.375" (6 cm)  
Overall Length: 7.375" (18,7 cm)

5470-20 [20 mm Chisel]

Osteotome Width: 20 mm  
Blade Length: 2.375" (6 cm)  
Overall Length: 7.375" (18,7 cm)

5472-08 [8 mm Offset  
Cement Removal Chisel]

Osteotome Dimensions: 8 mm Wide x 12 mm Long  
Blade Length: 2.375" (6 cm)  
Overall Length: 7.375" (18,7 cm)


5474-06 [6 mm Notched  
Cement Removal Chisel]

Osteotome Width: 6 mm  
Blade Length: 2.625" (6 cm)  
Overall Length: 7.375" (18,7 cm)

5475-08 [8 mm Implant Remover]

Diameter: 8 mm  
Blade Length: 2.625" (6 cm)  
Overall Length: 7.375" (18,7 cm)

5470-CASE [Case Only]



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.