


# Engh Unicondylar Minimally Invasive Knee Surgery Instruments

Designed by Gerard A. Engh, MD

*Designed for use in unicondylar minimally invasive knee surgery*



PRODUCT NO'S:	
4920-01	[Cement Scraper - Right] Scraper Head: 5 mm x 9 mm Overall Length: 8.5" (21,6 cm) Handle Length: 3.5" (8,9 cm)
4920-02	[Cement Scraper - Left] Scraper Head: 5 mm x 9 mm Overall Length: 8.5" (21,6 cm) Handle Length: 3.5" (8,9 cm)
4930-01	[Offset Osteotome] Osteotome Head: 10 mm x 10 mm Overall Length: 8.5" (21,6 cm) Handle Length: 3.5" (8,9 cm)

## *Cement Scrapers*

Right and left design used to scrape cement from around and behind the implants.

## *Offset Osteotome*

Designed to remove the overhanging posterior condyle or osteophytes from the back of the femoral condyle.

[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.