

# Evans Reverse Hohmann Retractor

Designed by Peter J. Evans, MD

*Designed for total shoulder arthroplasty and open rotator cuff procedures, the small size is useful for retracting the deltoid superiorly or laterally, and also protecting the axillary nerve inferiorly while simultaneously exposing the glenoid*



**PRODUCT NO:**

4547

Overall Length: 8.5" (21,6 cm)

Blade Width: Tapers from 30 mm to 18 mm

Blade Depth: 3" (7,6 cm)

Prong Width: 6 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**