

# Extra Deep Hip Retractors

All Extra Deep retractors are 2" (5 cm) longer than their standard version.

For hip surgery with large patients, and when extra large instruments are desired for increase depth and leverage



PRODUCT NO'S:
3418 [Extra Deep Mueller-type Femoral Neck Elevator] Overall Length: 15.25" (38,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width at Widest: 25 mm
4535-01 [Extra Deep Modified Narrow Hohmann Retractor] Overall Length: 11.5" (29,2 cm) Blade Width: 16.4 mm
4540-01 [Extra Deep Long Narrow Blunt Hohmann Retractor] Overall Length: 13.25" (33,7 cm) Blade Width: 22 mm Blade Width at End: 16 mm
4550-01 [Extra Deep Modified Blunt Hohmann Retractor] Overall Length: 13.25" (33,7 cm) Blade Width at End: 11 mm
4558-01 [Extra Deep Hohmann Retractor] Overall Length: 11.5" (29,2 cm) Blade Width: 16.7 mm
6450-01 [Extra Deep Single Prong Soft Tissue Retractor] Overall Length: 13.75" (34,9 cm) Blade Width: 22.3 mm
6450-04 [Extra Deep Single Prong Soft Tissue Retractor with Short Tip] Overall Length: 13" (33 cm) Blade Width: 22.3 mm
6570-01 [Extra Deep Single Prong Acetabular Retractor] Overall Length: 13.75" (34,9 cm) Blade Width: 22.3 mm
6595-01 [Extra Deep Modified Wide Hohmann Retractor] Overall Length: 11.5" (29,2 cm) Blade Width: 42.5 mm
7115-03 [Extra Deep Bent Hohmann Retractor] Overall Length: 12.125" (30,8 cm) Handle Length: 9.75" (24,8 cm) Depth from Bend: 6.25" (15,9 cm) Blade Width: 19 mm
7630-03 [Extra Deep Large Cobra Retractor] Overall Length: 19" (48,2 cm) Handle Length: 14" (35,6 cm) Blade Width at Widest: 33 mm

Extra Deep Mueller-type Femoral Neck Elevator design modification by Tom Eickmann, MD

- 1 Extra Deep Mueller-type Femoral Neck Elevator #3418
- 2 Extra Deep Modified Hohmann #4535-01
- 3 Extra Deep Long Narrow Blunt Hohmann #4540-01
- 4 Extra Deep Modified Blunt Hohmann #4550-01
- 5 Extra Deep Hohmann #4558-01
- 6 Extra Deep Single Prong Soft Tissue #6450-01
- 6S Extra Deep Single Prong Soft Tissue with Short Tip #6450-04
- 7 Extra Deep Single Prong Acetabular #6570-01
- 8 Extra Deep Modified Wide Hohmann #6595-01
- 9 Extra Deep Bent Hohmann #7115-03
- 10 Extra Deep Large Cobra #7630-03



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.