




Faillace Bone Impact/Graft Forceps

Design modification by John J. Faillace, MD, FAAOS

Long vertical grooves at the tip are designed to deliver graft into a small space, where a freer elevator can be used to push the graft down into the space, then the closed flat end can be used to tamp down the graft



PRODUCT NO:

5011

Overall Length: 5" (12,7 cm)

Tip Diameter When Closed: 3,2 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED