

Fajardo Femoral Neck Reference Retractor for Anterior THA

Designed by Jesse Fajardo, MD

Designed to be placed along the inferior femoral neck during anterior total hip arthroplasty to help expose the femoral neck, the retractor also has cutouts that allow the surgeon to mark the femoral neck for osteotomy at the desired level referencing the lesser trochanter, with the cutouts indicating 4 mm, 8 mm, 12 mm, and 16 mm

Can also be used in upper or lower extremity procedures where retracting and making measurements is needed, especially if fluoroscopy is used, because the known dimensions of the retractor and cutouts help provide scale to the fluoroscopic image.

2 mm cutouts every 4 mm

PRODUCT NO:

6419

Overall Length: 13" (33 cm)
Handle Length: 8.5" (21.6 cm)
Blade Width: .8" (2 cm)
Blade Drop: 5.5" (14 cm)
Blade Outcropping: 1" (2.54 cm)



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362