

# Femoral Tibial Coated Spreader Bar

*Designed to separate the femur and tibia when implant components are in place*

The end is coated to help protect from scratching component surfaces.



**PRODUCT NO:**

1820

Overall Length: 13" (33 cm)

Coated Surface: 4" (10,2 cm)

Blade Width: 13 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**