


# Femoral Head Removers

Used to remove a femoral head during total hip arthroplasty or fracture surgery



Quick-connect version for use with a driver.

#### PRODUCT NO'S:

3688 [Femoral Head Remover with Zimmer Hall Quick-connect] Overall Length: 8.5" (21,6 cm)
3690 [Femoral Head Remover with T-Handle] Overall Length: 8.75" (22,2 cm)
8248 [Fixed Driver with Zimmer Hall Quick-connect] Overall Length: 5.75" (15,6 cm) Handle Width: 4.625" (11,6 cm)



Fixed Driver  
with Zimmer Hall  
Quick-connect



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.